

Recommended cleaning and disinfection agents for upholstered surfaces made by Meden-Inmed

Akzo Nobel

Natriumhypochlorid

Bode

Bacillol AF

Cutasept

Braun

Meliseptol Foam Pure

Softasept N

Meliseptol

Helipur

Melsitt

Melsept

Desomed

Desomed Rapid

Dr. Schumacher

Descosept AF

Ecolab

Incidin Foam

Incidin Plus

Incidin Perfekt

Incidin Liquid

Favodent Karl Huber

Flora quick

Henkel

Incidur-Spray

Lysoform

Fugaten Spray

Orochemie

B45

Peppler

Micro clan

Pepp clean

Ruck

Sprühdesinfektion

Reckitt Benckiser

Sagrotan med

Schülke&Mayr

Kodan Tinktur Forte

Esemfix

Octenisept

Mikrozid AF

Unident

Unisepta foam

Please apply any of the agents first to a concealed area to test its influence on the upholstery.