

Meden-Inmed Sp. z o. o.
ul. Wenedów 2, 75-847 Koszalin, Poland


Phone.: +48 94 344 90 50, +48 94 347 10 52

Fax: +48 94 347 10 41

Mobile: +48 602 668 002, +48 601 646 866

E-mail: international@meden.com.pl

www.en.meden.com.pl



PREMIUM UPHOLSTERY



Management
System
ISO 13485:2016

www.tuv.com
ID 0000073249

2022.02.25/ WMD_4/ Wersja 4

Technical specification:

Coating: 100% Polyurethane
Phthalate - free

Treated with the most advanced silver ion technology, a natural environmentally friendly bacteriafighting defense to inhibit odor and stain causing bacteria.



6 times more durable than standard



Impermeable



Stretching resistant



Low temp. resistant



Extremely soft



UV resistant



Mildew resistant



Chemical resistant



Easy to clean



Abrasion resistant

Cleaning:

We recommend to clean this article with water and soap only.



6701 WHITE



6703 ECRU



6731 MANGO



6722 GRASS



6724 JET



6734 QUETZAL



6718 TANGELO




6738 WINE




6726 ARTESIAN



6739 GRAPHITE



6728 BLACK



6713 CHOCOLATE