




# Massage table catalogue

## 2023/2024

The catalogue contains information on, among other things, medical devices, that should be operated by qualified medical personnel and used in accordance with the user manual or label.

#### Remarks:

- Medical devices are marked with sign 
- In addition to massage tables, you will also find here other types of tables, e.g. tilting or examination and treatment tables.
- This catalogue is for information purposes and cannot be treated as a commercial offer within the meaning of Art. 66 CC.
- Manufacturer hereby reserves the right to change design that would not disturb the principal functionality and safety requirements. Despite every effort, the colours presented in the catalogue can be differ from the real due to colour printing rendition reasons.
- Images in this catalogue may contain options that are subject to a fee.



Meden-Inmed was founded in 1989 as a company employing two people. Since then, it has gained a strong market position through hard work and reliability. Initially a small distributor of a few renowned manufacturers, now a dynamic company with over 470 employees - professionals in the fields of electronics, mechanics, economy, sales and marketing.

Modern rehabilitation combines specialized knowledge with state of the art technology. We knew that when we started our own production plant in 1996 and are even more convinced of it now, when we create new products and widen our offer of medical devices each year. We make sure they are modern, efficient, comfortable and user friendly. We are always happy to see new ideas turn into new products.

Every product we have sold proves the trust our customers have placed in us. Behind every device we deliver there is a long-term relationship. Our customers consciously select a supplier that not only delivers the equipment, but can also guarantee quality and professional service. It has always been our duty and our motto.

To address the growing needs of our customers we have now expanded the range of our products. Our goal is to deliver you complete range of equipment required by professionals in the field of rehabilitation. This catalogue shows just a part of our offer. We encourage you to visit our website.

We will help you choose the best equipment that is tailored to your needs. Let this catalogue be your guide to our products. Meden-Inmed's Team will assist you in every matter of our cooperation. We are always here to hear your voice. Your feedback is always welcome and appreciated.

**Maciej Zinka**  
CEO

**Ph.D. Eng. Wiesław Zinka**  
Senior President

## Our awards:



**European Medal 2023, 2022, 2021, 2020, 2016, 2013, 2012, 2011, 2010, 2009, 2008, 2007**

European Economic and Social Committee and Business Center Club



**Diamond of FORBES Magazine 2009, 2019**



**Ambassador of Polish Economy in category:**

**Exporter, 2014, 2022  
XXI century solution creator, 2013  
Partner of foreign companies, 2012**



**Perfect Line 2009 SPA Business**



**Gazelle of Business 2008  
Gazelle of Business 2007  
Gazelle of Business 2006**



**Awards of biggest rehabilitation equipment fairs in Poland - Rehabilitacja Fairs**

**XXVII Rehabilitation Fair Łódź 2019  
XXII Rehabilitation Fair Łódź 2014  
XXI Rehabilitacja Łódź 2013  
XVI Rehabilitacja Łódź 2008  
XV Rehabilitacja Łódź 2007  
XIV Rehabilitacja Łódź 2006  
XIII Rehabilitacja Łódź 2005  
X Rehabilitacja Łódź 2002**



**5th rank on the list of 100 best companies in the region**

II Euroregional Economic Presentation Szczecin 2002



**Investment of the year 2009**

Polish Federation of Engineering Associations



**Product of West Pomerania 2015, 2012, 2010, 2009  
Product of Pomerania 2009, 2008, 2007**

Polish Federation of Engineering Associations



**Product of Koszalin 2014, 2013, 2010, 2009, 2008, 2007, 2006**

Polish Federation of Engineering Associations



**Best Distributor Award 2019**

KOELIS



**Pioneer Distributor Award 2017**

Business Centre Club




**Leader of Polish business 2016, 2017, 2019**

Business Center Club

## Types of wheeling system:


### F0

Composed of 2 legs with rubber feet and 2 small wheels for moving of Prestige Table by lifting it from one side. It is recommended for users who relocate their Prestige Table less frequently.



## Massage tables


### Terapeuta SCM Prestige





## F4

Composed of 4 directional castors with central brake allows you very easy and fast relocation of Prestige Table with or without a patient.



Two-layer padding with upholstery available in many colours



Release handle for gas spring system in leg section



Foot - operated lever of F4 system



Rubber foot with 1 cm height adjustment



# Years of experience

..., our individual research and passion for innovations led to creation of a stable, **reliable and extremely functional table line** for massage.

You can now enjoy the **high quality and high functionality of Therapie Prestige Table Line**. Just choose one of Prestige Tables which will suit your needs and the needs of your patients. We care about your opinion so if you think that something should be changed or added in the current Therapie Prestige Table Line please let us know what you think.

## Unique features of Prestige Tables:

### 1. High quality upholstery!

In order to guarantee comfort of patients and high functionality for users each table has a **two-layer padding** with specially selected density of foams. After years of checking various upholstery samples, we finally found an upholstery which is:

- Non-flammable
- Bio-compatible
- Scratchproof



High quality finish

The features of upholstery together with high durability of the powder paint used for covering the metal elements will allow you to enjoy your Therapie Prestige Table in your practice.


There is possibility to buy couches also with upholstery made with 100% Polyurethane (**Premium Upholstery** or **OMEGA Upholstery**).

### 2. We care about your working comfort!

That's why we invented the unique **HalloTronic®** System that enables high/low height adjustment via hallotrons without the use of mechanical switches. Thanks to **HalloTronic®** System you can choose the right height of your Prestige Table without effort and noise. This solution assures the comfort of patient and comfort of work. **HalloTronic®** System allows you to adjust height of your table with help of a foot bar from each side of your Prestige Table with just a slight movement of your foot. Enjoy the benefits of your Prestige Table and unique features of **HalloTronic®** System:

- High/low height adjustment via electronically controlled actuator
- No noise during lifting
- Long life span
- Reliability
- Operated with minimal force by foot
- High/low height adjustment via shifting of bar, easily accessible from each side of Prestige Table

Prestige Tables ensure smooth high/low height adjustment with a powerful electric actuator with lifting capacity up to 225 kg.



### 3. Wheeling systems tailored to your needs!

You can choose between F0 and F4 wheeling system. F0 system composed of 2 legs and 2 small wheels allows for moving of Prestige Table by lifting it from one side. The system is recommended for users who relocate their Prestige Table less frequently. F4 system composed of 4 directional castors with central brake allows you for a very easy and fast relocation of Prestige Table with or without a patient.

### 4. Easy adjustment as you wish!


Enjoy the **easy adjustment of Prestige Table's sections** with the help of gas springs. Ergonomically located and designed release handles **allow you to adjust sections of your Prestige Table the way you need** practically without any effort and with one hand. Choose postural drainage (pivot) option to allow for easy positioning of the patient in order to perform better. You are able to adjust the postural drainage via an electronically controlled actuator by easily accessible two postural drainage control switches from both sides of your Prestige Table. Unique self-locking fold-down armrest (SmartS4Arms) for a very easy (one hand operated) adjustment.



### 5. Guarantee your practice the ultimate safety!

All Prestige Tables have an integrated electronic control unit allowing for height adjustment and adjustment of the postural drainage option with an indication of stand-by mode (green light) and **Personal Authorization Safety System** (PASS). Each Prestige Table is equipped with two magnetic safety keys (PASS Keys) in order to prevent unauthorized use of a table which may lead to severe accidents and compensation claims. What is more, all Prestige Tables have a **special safety system** that prevents table from folding down in case of power shortage or actuator damage.

Tables are manufactured in conformity with the certified quality management system: ISO 13485:2016.



### 6. The devil is in the details!

All details of Prestige Tables are designed to help you in your practice, with their functionality, durability and great design. You can take advantage of **semi-automatic blocking/release system for fold-down side sections, comfortable armrests and ergonomic table top** allowing for an easy access to the patient. **Rubber feet** increase the stability. The patients can feel comfort of the **anatomical design of the face hole**. We covered screws and other protruding elements with unobtrusive plastic covers.

### 7. We care about the production process!


All frames are welded automatically by high quality robotic system, which makes Prestige Tables **very stable and resistant over years**. Metal parts are being covered using the most advanced technological processes in a modern powder painting line. This means that you will enjoy the high quality of your Prestige Table for years and the maintenance (e.g. cleaning) will be a pleasure.



# Terapeuta SCM-2 Prestige M-S2.F0/M-S2.F4

Massage table

MD



Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**



Length: **202 cm**



Height: **47-95 cm**



Width: **69 cm**



Weight: **80/95 kg**



Max load\*: **225 kg**

## Features:

- 2-section massage table.
- Adjustable head section from -70° to +40° with a gas spring.
- Face hole in the head section with 2 removable face hole covers (round, crescent).
- Two-layer padding with upholstery in many colours.
- Powder coated, extremely stable steel frame.
- 2 side brackets for fixation of stabilization belts and elastic bands.
- Traction support sleeves.
- Personal Authorization Safety System with 2 safety magnetic keys.
- HalloTronic® System for electric high/low height adjustment from 47 to 95 cm with a lifting capacity up to 225 kg.
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table.
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.
- Two ways of frame painting: „white” or „twintone”.


## Additional option:

- Built-in adjustable surface heating system.
- Upholstery of 80 mm thickness.
- Upholstery of 80 cm width.

## Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Traction supports.
- Stabilization belts.
- Stools for therapist.

## Table articulations:



## Ordering information:

**M-S2.F0:** 2-section massage table with wheeling system composed of 2 small wheels and 2 legs with rubber feet.

**M-S2.F4:** 2-section massage table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.

\* max. lifting capacity



# Terapeuta SCM-3 Prestige M-S3.F0/M-S3.F4

Massage table

MD



Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**



Length: **205 cm**



Height: **47-95 cm**



Width: **69 cm**



Weight: **85/100 kg**



Max load\*: **225 kg**

## Features:

- 3-section massage table.
- Adjustable head section from -70° to +40° with a gas spring.
- Face hole in the head section with 2 removable face hole covers (round, crescent).
- Adjustable leg section from -10° to 80° with 2 gas springs.
- Release handle of the gas spring system in the leg section accessible from both sides of the table.
- Two-layer padding with upholstery in many colours.
- Powder coated, extremely stable steel frame.
- Traction support sleeves.
- Personal Authorization Safety System with 2 safety magnetic keys.
- HalloTronic® System for electric high/low height adjustment from 47 to 95 cm with a lifting capacity up to 225 kg.
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table.
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.
- Two ways of frame painting: „white” or „twintone”.

## Additional option:

- Built-in adjustable surface heating system.
- Upholstery of 80 mm thickness.
- Upholstery of 80 cm width.

## Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Traction supports.
- Stools for therapist.


## Ordering information:

**M-S3.F0:** 3-section massage table with wheeling system composed of 2 small wheels and 2 legs with rubber feet.

**M-S3.F4:** 3-section massage table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.

\* max. lifting capacity

## Table articulations:






# Terapeuta SCM-3

## Prestige M-P3.F0/M-P3.F4

Massage table

MD



2 ways of frame painting



Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**



Length: **205 cm**



Height: **51-99 cm**



Width: **69 cm**



Weight: **95/110 kg**



Max load\*: **225 kg**

### Features:

- 3-section massage table.
- Postural drainage with an electric actuator operated by 2 switches from both sides of the table.
- Adjustable head section from -70° to +40° with a gas spring.
- Face hole in the head section with 2 removable face hole covers (round, crescent).
- Adjustable leg section from -10° to 80° with a gas spring.
- Release handle of the gas spring system in the leg section accessible from both sides of the table.
- Two-layer padding with upholstery in many colours.
- Powder coated, extremely stable steel frame.
- Traction support sleeves.
- Personal Authorization Safety System with 2 magnetic safety keys.
- HalloTronic® System for electric high/low height adjustment from 51 to 99 cm with a lifting capacity up to 225 kg.
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table.
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.
- Two ways of frame painting: „white” or „twintone”.


### Additional option:

- Built-in adjustable surface heating system.
- Upholstery of 80 mm thickness.
- Upholstery of 80 cm width.

### Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Traction supports.
- Stools for therapist.

### Table articulations:



### Ordering information:

**M-P3.F0:** 3-section massage table with wheeling system composed of 2 small wheels and 2 legs with rubber feet.

**M-P3.F4:** 3-section massage table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.

\* max. lifting capacity




# Terapeuta SCM-4

## Prestige M-S4.F0/M-S4.F4

Massage table



MD



2 ways of frame painting



Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**



Length: **202 cm**



Height: **47-95 cm**



Width: **69 cm**



Weight: **87/102 kg**



Max load\*: **225 kg**

### Features:

- 4-section massage table.
- Adjustable head section from -70° to +40° with a gas spring.
- Face hole in the head section with 2 removable face hole covers (round, crescent).
- Innovative self-locking fold-down armrests in head section (one hand operated) - SmartS4Arms (Smart Slide For Arms).
- Two-layer padding with upholstery in many colours.
- Powder coated, extremely stable steel frame.
- 2 side brackets for the fixation of stabilization belts and elastic bands.
- Traction support sleeves.
- Personal Authorization Safety System with 2 magnetic safety keys.
- HalloTronic® System for electric high/low height adjustment from 47 to 95 cm with a lifting capacity up to 225 kg.
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table.
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.
- Two ways of frame painting: „white” or „twintone”.

### Additional option:

- Built-in adjustable surface heating system.
- Upholstery of 80 mm thickness.
- Upholstery of 80 cm width.

### Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Traction supports.
- Stabilization belts.
- Stools for therapist.


### Ordering information:

**M-S4.F0:** 4-section massage table with wheeling system composed of 2 small wheels and 2 legs with rubber feet.

**M-S4.F4:** 4-section massage table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.

\* max. lifting capacity

### Table articulations:



# Terapeuta SCM-5 Prestige M-S5.F0/M-S5.F4

Massage table



MD



2 ways of frame painting



Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**



Length: **205 cm**



Height: **47-95 cm**



Width: **69 cm**



Weight: **92/107 kg**



Max load\*: **225 kg**

## Features:

- 5-section massage table.
- Adjustable head section from -70° to +40° with a gas spring.
- Face hole in the head section with 2 removable face hole covers (round, crescent).
- Innovative self-locking fold-down armrests in head section (one hand operated) - SmartS4Arms (Smart Slide For Arms).
- Adjustable leg section from -10° to 80° with 2 gas springs.
- Release handle of the gas spring system in the leg section accessible from both sides of the table.
- Two-layer padding with upholstery in many colours.
- Powder coated, extremely stable steel frame.
- Personal Authorization Safety System with 2 magnetic safety keys.
- HalloTronic® System for electric high/low height adjustment from 47 to 95 cm with a lifting capacity up to 225 kg.
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table.
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.
- Two ways of frame painting: „white” or „twintone”.


## Additional option:

- Built-in adjustable surface heating system.
- Upholstery of 80 mm thickness.
- Upholstery of 80 cm width.

## Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Traction supports.
- Stools for therapist.

## Table articulations:



## Ordering information:

**M-S5.F0:** 5-section massage table with wheeling system composed of 2 small wheels and 2 legs with rubber feet.

**M-S5.F4:** 5-section massage table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.

\* max. lifting capacity




# Terapeuta SCM-5

## Prestige M-P5.F0/M-P5.F4

Massage table



MD



2 ways of frame painting



Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**



Length: **205 cm**



Height: **51-99 cm**



Width: **69 cm**



Weight: **102/117 kg**



Max load\*: **225 kg**

### Features:

- 5-section massage table.
- Postural drainage with an electric actuator operated by 2 switches from both sides of the table.
- Adjustable head section from -70° to +40° with a gas spring.
- Face hole in the head section with 2 removable face hole covers (round, crescent).
- Innovative self-locking fold-down armrests in head section (one hand operated) - SmartS4Arms (Smart Slide For Arms).
- Adjustable leg section from -10° to 80° with a gas spring.
- Release handle of the gas spring system in the leg section accessible from both sides of the table.
- Two-layer padding with upholstery in many colours.
- Powder coated, extremely stable steel frame.
- Personal Authorization Safety System with 2 magnetic safety keys.
- HalloTronic® System for electric high/low height adjustment from 51 to 99 cm with a lifting capacity up to 225 kg.
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table.
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.
- Two ways of frame painting: „white” or „twintone”.


### Additional option:

- Built-in adjustable surface heating system.
- Upholstery of 80 mm thickness.
- Upholstery of 80 cm width.

### Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Traction supports.
- Stools for therapist.

### Table articulations:



### Ordering information:

**M-P5.F0:** 5-section massage table with wheeling system composed of 2 small wheels and 2 legs with rubber feet.

**M-P5.F4:** 5-section massage table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.

\* max. lifting capacity




# Terapeuta SCM-7

## Prestige M-S7.F0/M-S7.F4

Massage table



MD



2 ways of frame painting



Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**



Length: **205 cm**



Height: **47-95 cm**



Width: **69 cm**



Weight: **95/110 kg**



Max load\*: **225 kg**

### Features:

- 7-section for massage table.
- Adjustable head section from -70° to +40° with a gas spring.
- Face hole in the head section with 2 removable face hole covers (round, crescent).
- Fold-down side sections in the middle section with a semi-automatic blocking/release system.
- Innovative self-locking fold-down armrests in head section (one hand operated) - SmartS4Arms (Smart Slide For Arms).
- Adjustable leg section from -10° to 80° with a gas spring.
- Release handle of the gas spring system in the leg section accessible from both sides of the table.
- Two-layer padding with upholstery in many colours.
- Powder coated, extremely stable steel frame.
- Traction support sleeves.
- Personal Authorization Safety System with 2 magnetic safety keys.
- HalloTronic® System for electric high/low height adjustment from 47 to 95 cm with a lifting capacity up to 225 kg.
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table.
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.
- Two ways of frame painting: „white” or „tintone”.


### Additional option:

- Built-in adjustable surface heating system.
- Upholstery of 80 mm thickness.
- Upholstery of 80 cm width.

### Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Traction supports.
- Stools for therapist.

### Table articulations:



### Ordering information:

**M-S7.F0:** 7-section massage table with wheeling system composed of 2 small wheels and 2 legs with rubber feet.

**M-S7.F4:** 7-section massage table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.

\* max. lifting capacity






# Terapeuta SCM-7

## Prestige M-P7.F0/M-P7.F4

Massage table



MD



2 ways of frame painting



Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**



Length: **205 cm**



Height: **51-99 cm**



Width: **69 cm**



Weight: **105/120 kg**




Max load\*: **225 kg**

### Features:

- 7-section massage table.
- Postural drainage with an electric actuator operated by 2 switches from both sides of the table.
- Adjustable head section from -70° to +40° with a gas spring.
- Face hole in the head section with 2 removable face hole covers (round, crescent).
- Innovative self-locking fold-down armrests in head section (one hand operated) - SmartS4Arms (Smart Slide For Arms).
- Adjustable leg section from -10° to 80° with a gas spring.
- Fold-down side sections in the middle section with a semi-automatic blocking/release system.
- Release handle of the gas spring system in the leg section accessible from both sides of the table.
- Two-layer padding with upholstery in many colours.
- Powder coated, extremely stable steel frame.
- Traction support sleeves.
- Personal Authorization Safety System with 2 magnetic safety keys.
- HalloTronic® System for electric high/low height adjustment from 51 to 99 cm with a lifting capacity up to 225 kg.
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table.
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.
- Two ways of frame painting: „white” or „twintone”.

### Table articulations:



### Additional option:

- Built-in adjustable surface heating system.
- Upholstery of 80 mm thickness.
- Upholstery of 80 cm width.

### Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Traction supports.
- Stools for therapist.

### Ordering information:

**M-P7.F0:** 7-section massage table with wheeling system composed of 2 small wheels and 2 legs with rubber feet.

**M-P7.F4:** 7-section massage table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.


\* max. lifting capacity



# Terapeuta SCM-2 Prestige E-S2.F0/E-S2.F4

Massage table

MD



2 ways of frame painting

Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**

Length: **202 cm**

Height: **47-95 cm**

Width: **69 cm**

Weight: **80/95 kg**

Max load\*: **225 kg**

## Features:

- 2-section massage table with a long head section.
- Adjustable long head section from 0° to +85° with 2 gas springs.
- Face hole in the head section with a removable face hole cover.
- Two-layer padding with upholstery in many colours.
- Powder coated, extremely stable steel frame.
- Optional 2 side brackets for the fixation of stabilization belts and elastic bands or 2 retractable side railings for fall prevention.
- Personal Authorization Safety System with 2 safety magnetic keys.
- HalloTronic® System for electric high/low height adjustment from 47 to 95 cm with a lifting capacity up to 225 kg
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table.
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.
- Two ways of frame painting: „white” or „twintone”.


## Additional option:

- Built-in adjustable surface heating system.
- Upholstery of 80 mm thickness.
- Upholstery of 80 cm width.

## Additional accessories:

- Paper holder.
- Optional: 2 side brackets or 2 retractable side railings.
- Wedges.
- Bolsters, semi-bolsters.
- Stabilization belts.
- Stools for therapist.

## Table articulations:



## Ordering information:

**E-S2.F0:** 2-section massage table with wheeling system composed of 2 small wheels and 2 legs with rubber feet.

**E-S2.F4:** 2-section massage table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.

\* max. lifting capacity




# Terapeuta SCM-5

## Prestige O-S5.F0/O-S5.F4

Massage table

MD



2 ways of frame painting



Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**



Length: **205 cm**



Height: **47-95 cm**



Width: **59 cm**



Weight: **92/107 kg**



Max load\*: **225 kg**

### Features:

- 5-section table.
- Narrow (31 cm), adjustable head section from -70° to +40° with gas spring.
- Face hole in the head section with a removable face hole cover.
- Fold-down side sections allowing broadening of middle section if necessary.
- Adjustable leg section from -10° to +80° with 2 gas springs.
- Release handle of the gas spring system in the leg section accessible from both sides of the table.
- Two-layer padding with upholstery in many colours.
- Powder coated, extremely stable steel frame.
- Traction support sleeves.
- Personal Authorization Safety System with 2 magnetic safety keys.
- HalloTronic® System for electric high/low height adjustment from 47 to 95 cm with a lifting capacity up to 225 kg.
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table.
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.
- Ergonomic design.
- Two ways of frame painting: „white” or „twintone”.


### Additional option:

- Built-in adjustable surface heating system.
- Upholstery of 80 mm thickness.

### Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Traction supports.
- Stools for therapist.

### Table articulations:



### Ordering information:

**O-S5.F0:** 5-section massage table with wheeling system composed of 2 small wheels and 2 legs with rubber feet.

**O-S5.F4:** 5-section massage table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.

\* max. lifting capacity



# Terapeuta SCM-5 Prestige O-P5.F0/O-P5.F4

Massage table

MD



2 ways of frame painting



Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**



Length: **205 cm**



Height: **51-99 cm**



Width: **59 cm**



Weight: **102/117 kg**



Max load\*: **225 kg**

## Features:

- 5-section table.
- Postural drainage with electric actuator operated by 2 switches from both sides of the table.
- Narrow (31 cm), adjustable head section from -70° to +40° with gas spring.
- Face hole in the head section with a removable face hole cover.
- Fold-down side sections allowing for the broadening of middle section if necessary.
- Adjustable leg section from -10° to +80° with 2 gas springs.
- Release handle of the gas spring system in the leg section accessible from both sides of the table.
- Two-layer padding with upholstery in many colours.
- Powder coated, extremely stable steel frame.
- Traction support sleeves.
- Personal Authorization Safety System with 2 magnetic safety keys.
- HalloTronic® System for electric high/low height adjustment from 51 to 99 cm with a lifting capacity up to 225 kg.
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table.
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.
- Ergonomic design.
- Two ways of frame painting: „white” or „twintone”.


## Additional option:

- Built-in adjustable surface heating system.
- Upholstery of 80 mm thickness.

## Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Traction supports.
- Stools for therapist.

## Table articulations:



## Ordering information:

**O-P5.F0:** 5-section table with wheeling system composed of 2 small wheels and 2 legs with rubber feet.

**O-P5.F4:** 5-section table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.

\* max. lifting capacity



# Terapeuta Prestige Line




Technical data:											
	M-S2	M-S3	M-P3	M-S4	M-S5	M-P5	M-S7	M-P7	E-S2	O-S5	O-P5
Dimensions (length x width) [cm]:	202 x 69	205 x 69	205 x 69	202 x 69	205 x 69	205 x 69	205 x 69	205 x 69	202 x 69	205 x 59	205 x 59
Number of sections:	2	3	3	4	5	5	7	7	2	5	5
Head section inclination angle [°]:	-70/+40	-70/+40	-70/+40	-70/+40	-70/+40	-70/+40	-70/+40	-70/+40	0/+85	-70/+40	-70/+40
Middle section inclination angle [°]:	-	-	0/+24	--	-	0/+24	-	0/+24	-	-	0/+24
Leg section inclination angle [°]:	-	-10/+80	-10/+80	-	-10/+80	-10/+80	-10/+80	-10/+80	-	-10/+80	-10/+80
Height adjustment [cm]:	47-95	47-95	51-99	47-95	47-95	51-99	47-95	51-99	47-95	47-95	51-99
Weight F0/F4 [kg]:	80/95	85/100	95/110	87/102	92/107	102/117	95/110	105/120	80/95	92/107	102/117
Max. lifting capacity [kg]:	225	225	225	225	225	225	225	225	225	225	225
Power supply [V/Hz]:	100-240/50-60	100-240/50-60	100-240/50-60	100-240/50-60	100-240/50-60	100-240/50-60	100-240/50-60	100-240/50-60	100-240/50-60	100-240/50-60	100-240/50-60
Max. working load [kg]:	300	300	300	300	300	300	300	300	300	300	300



Powerful electric actuator



Adjustment of postural drainage



Breathing slot shaped to place the nose and chin



Personal Authorization Safety System (PASS)




Hallotronic® System for electric height adjustment



## Types of base:

### F0

Composed of 2 legs with rubber feet and 2 small wheels for moving of Elemental Table by lifting it from one side. The system is recommended for users who relocate their Elemental Table less frequently.




## Massage tables Elemental



## F0 + F2

Made on two small wheels and two feet, with an attached mechanism F2 of two swivel wheels, facilitating table transport.



Hand controller

Two-layer padding with upholstery in many colours



Transport wheels


Powerful electric actuator



Electric height adjustment



SmartS4Arms



High/low height adjustment foot bar - option  
(replaces foot- or hand-operated controller)



# Elemental Ice

Massage table


## Tables articulations:

ICE S2.F0:

MD



(Hand-operated controller)





ICE S4.F0:

MD



(Hand-operated controller)





ICE E2.F0:

MD



(Foot-operated controller)




# Elemental Fire

Massage table

## Tables articulations:


### FIRE S3.F0:

MD



### FIRE S5.F0:

MD



## Additional accessories:


### Lifter frame F2




### Horizontal folding barriers (E2 only)




### Folding barriers (E2 only)



### Paper holder type 3



### Foot- or hand-operated controller



### Wedges, bolsters, traction supports, belts



# Elemental Ice, Fire

Massage table

MD



Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**



Length: **202-205 cm**



Height: **48-96 cm**



Width: **69 cm**



Weight: **80-88 kg**



Max load\*: **180 kg**

## Features:

- 2-, 3-, 4-, 5-section massage tables.
- Adjustable head, leg section with gas springs.
- Round face hole in the head section with 2 removable face hole covers (crescent, round).
- Self-locking fold-down armrest in head section (S4, S5) - SmartS4Arms (Smart Slide For Arms).
- Two-layer padding with upholstery in many colours.
- Powder coated, extremely stable two-colour steel frame.
- 2 side brackets for the fixation of stabilization belts.
- Electric height adjustment: 48-96 cm via hand-operated remote controller (with authorization) or foot controller to choose in standard.
- High/low height adjustment foot bar accessible from each side of the table - option (replaces remote controller).
- Lifting capacity up to 180 kg, working load up to 250 kg.
- Locking safety system on the remote controller.
- Basic wheeling system F0 composed of 2 small wheels and 2 legs with rubber feet.
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.



Crescent cover

## Technical data:

	ICE S2	ICE E-S2	ICE S4	FIRE S3	FIRE S5
Dimensions (length x width) [cm]:	202 x 69	202 x 69	202 x 69	205 x 69	205 x 69
Number of sections:	2	2	4	3	5
Head section inclination angle [°]:	-70/+40	-20/+85	-70/+40	-70/+40	-70/+40
Middle section inclination angle [°]:	--	--	--	--	--
Leg section inclination angle [°]:	--	--	--	0/+85	0/+85
Height adjustment [cm]:	48-96	48-96	48-96	48-96	48-96
Weight [kg]:	81,5	81	80	86	88
Max. lifting capacity [kg]:	180	180	180	180	180
Power supply [V/Hz]:	100-240/50-60	100-240/50-60	100-240/50-60	100-240/50-60	100-240/50-60
Max. working load [kg]:	250	250	250	250	250

\* max. lifting capacity





# Accessories



Pillow



Traction support



Stabilization belt



Kaltenborn wedge (small/large)



Wedges (different sizes)



Semi-bolsters (different sizes)



Bolsters (different sizes)



Magnetic safety key




Covers for face hole

## Upholstery colour options:

Enjoy the high quality upholstery in the colours you like!

You can now choose non-flammable, bio-compatible and scratchproof upholstery for your tables, couches, stools as well as wedges and bolsters in colours matching the design and atmosphere of your place. Softness and durability of the upholstery provide comfort to patients as well as therapists.

## Standard upholstery colour:



## Paper holders for:

**PRESTIGE E-S2, SAFARI E-S2  
THERASPACE E-S2  
ELEMENTAL E2**

type 1  
to upper frame



R9002DS  
ANTRACYT

**ELEMENTAL**

type 3  
to lower frame



R9007DS

**SAFARI, PRESTIGE**

type 8  
to lower frame



R9006DS

**SAFARI ELEPHANT**

type 8B  
to lower frame



R9006DS

**NEXUS**

type 12  
to upper frame



R9002DS

**THERASPACE M-S2, M-S4**

type 4  
to upper frame



R9007DS

**THERASPACE X, Z**

type 13  
to upper frame



R9006DS

## Premium upholstery colour:

Topcoat: 100% Polyurethane



There is possibility to order (as option) also in other polyurethane upholstery - OMEGA - see p. 35



## Wheeling system:


### F4

Composed of 4 liftable directional wheels with a central locking mechanism and 4 rubber feet.



## Massage tables


### Theraspace





Sections adjust  
control switches

High-quality, soft and  
pleasant to the touch  
polyurethane upholstery



Release handle for  
gas spring system  
in leg section

F4 system central  
wheel locking  
mechanism lever

Base combining  
wood and steel

High/low height  
adjustment foot bar

Directional wheels facilitating  
the table transport



# Theraspace E-S2.F4

Massage table



MD



High/low height adjustment foot bar



Box dimensions: **220 x 85 x 80 cm** Weight: **+50 kg**



Length: **202 cm**



Height: **47-95 cm**



Width: **69 cm**



Weight: **115 kg**



Load\*: **200 kg**

## Features:

- 2-section massage table with long head section.
- High-quality, soft and pleasant to the touch polyurethane upholstery.
- Two-layer padding (2 cm foam 40 kg/cm<sup>3</sup> + 2 cm foam 44 kg/cm<sup>3</sup>).
- Stretched upholstery.
- Adjusting the angle of the head section with a gas spring (from -20 to 85°).
- Round breathing hole in head section with 2 types of plugs (crescent, round).
- Total length of the bed 202 cm.
- Maximum safe load of the table when raising-lowering the table ≤ 200 kg.
- Maximum safe working load of the table ≤ 275 kg.
- Safe electric actuator against crushing.
- High/low height adjustment foot bar (47-95 cm).
- F4 wheel system consisting of 4 liftable directional wheels with a central locking mechanism and 4 feet.
- A frame combining wood and steel.
- Protection against accidental start-up.
- Linak actuator.


## Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Horizontal folding barriers.
- Stools for therapist.



Stretched upholstery

## Settings and dimensions:



## Ordering information:

**E-S2.F4:** 2-section massage table with wheeling system composed of 4 liftable directional wheels with a central locking mechanism and 4 rubber feet.

\* max. safe load when raising-lowering the table



# Theraspace M-S2.F4 | Massage table



MD



High/low height adjustment foot bar



Box dimensions: 220 x 85 x 80 cm Weight: +50 kg



Length: 203 cm



Height: 47-95 cm



Width: 69 cm



Weight: 115 kg




Load\*: 200 kg

## Features:

- 2-section massage table.
- High-quality, soft and pleasant to the touch polyurethane upholstery.
- Two-layer padding (2 cm foam 40 kg/cm<sup>3</sup> + 2 cm foam 44 kg/cm<sup>3</sup>).
- Stretched upholstery.
- Round breathing hole in head section with 2 types of plugs (crescent, round).
- Adjusting the angle of the head section with a gas spring (from -70° to 40°).
- Total length of the bed 203 cm.
- Maximum safe load of the table when raising-lowering the table ≤ 200 kg.
- Maximum safe working load of the table ≤ 275 kg.
- Safe electric actuator against crushing.
- 2 side brackets for fixation of stabilization belts.
- 2 traction support sleeves.
- High/low height adjustment foot bar (47-95 cm).
- F4 wheel system consisting of 4 liftable directional wheels with a central locking mechanism and 4 rubber feet.
- A frame combining wood and steel.
- Protection against accidental start-up.
- Linak actuator.


## Additional accessories:

- Paper holder.
- Wedges.
- Traction supports.
- Stabilization belts.
- Bolsters, semi-bolsters.
- Stools for therapist.



Stretched upholstery

## Settings and dimensions:



## Ordering information:

**M-S2.F4:** 2-section massage table with wheeling system composed of 4 liftable directional wheels with a central locking mechanism and 4 rubber feet.

\* max. safe load when raising-lowering the table



# Theraspace M-S4.F4 | 2-section massage table



MD



High/low height adjustment foot bar



Box dimensions: **220 x 85 x 80 cm** Weight: **+50 kg**



Length: **203 cm**



Height: **47-95 cm**



Width: **69 cm**



Weight: **115 kg**




Load\*: **200 kg**

## Features:

- 4-section massage table.
- High-quality, soft and pleasant to the touch polyurethane upholstery.
- Two-layer padding (2 cm foam 40 kg/cm<sup>3</sup> + 2 cm foam 44 kg/cm<sup>3</sup>).
- Stretched upholstery.
- Round breathing hole in head section with 2 types of plugs (crescent, round).
- Adjusting the angle of the head section with a gas spring (from -70° to 40°).
- Total length of the bed 203 cm.
- Maximum safe load of the table when raising-lowering the table ≤ 200 kg.
- Maximum safe working load of the table ≤ 275 kg.
- Safe electric actuator against crushing.
- High/low height adjustment foot bar (47-95 cm).
- 2 side brackets for fixation of stabilization belts.
- 2 traction support sleeves.
- F4 wheel system consisting of 4 liftable directional wheels with a central locking mechanism and 4 rubber feet.
- A frame combining wood and steel.
- Protection against accidental start-up.
- Linak actuator.
- One-handed adjustment of armrest - **Smart4Arms**



## Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Traction supports.
- Stabilization belts.
- Stools for therapist.



Stretched upholstery

## Settings and dimensions:



## Ordering information:

**M-S4 .F4:** 4-section massage table with wheeling system composed of 4 liftable directional wheels with a central locking mechanism and 4 rubber feet.

\* max. safe load when raising-lowering the table





# Theraspace Z3.F4 | Massage table



MD



Box dimensions: **220 x 85 x 80 cm** Weight: **+50 kg**



Length: **205 cm**



Height: **49-95 cm**



Width: **67 cm**



Weight: **110 kg**



Load\*: **200 kg**

## Features:

- 3-section massage table.
- High-quality, soft and pleasant to the touch polyurethane upholstery.
- Progressive foam.
- Two-layer padding (3 cm foam 90 kg/cm<sup>3</sup> + 1 cm foam 43 kg/cm<sup>3</sup>).
- Stitched upholstery.
- Round breathing hole in head section with 2 types of plugs (crescent, round).
- Adjusting the angle of the head section with a gas spring (from -65 to 30°).
- Middle section adjustment angle: -20°.
- Head + middle section adjustment angle: +8°.
- Leg section adjustment angle: +8°.
- Total length of the table 205 cm.
- Maximum safe load of the table when raising-lowering the table ≤ 200 kg.
- Maximum safe working load of the table ≤ 275 kg.
- Safe electric actuator against crushing.
- High/low height adjustment foot bar (49-95 cm).
- F4 wheel system consisting of 4 liftable directional wheels with a central locking mechanism and 4 rubber feet.
- A frame combining wood and steel.
- Wireless remote control with three programmable positions (option).
- Electric adjustment of: head + middle, middle + leg sections via control switches.
- Protection against accidental start-up.
- Linak actuator.

## Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Stools for therapist.

## Additional option:


- Wireless remote control.

## Ordering information:

**Z3.F4:** 3-section massage table with wheeling system composed of 4 liftable directional wheels with a central locking mechanism and 4 rubber feet.

\* max. safe load when raising-lowering the table

## Settings and dimensions:



Remote control



Stitched upholstery



# Theraspace Z5.F4 | Massage table



MD



Box dimensions: **220 x 85 x 80 cm** Weight: **+50 kg**

Length: **205 cm**

Height: **49-95 cm**

Width: **67 cm**


Weight: **110 kg**

Load\*: **200 kg**

## Features:

- 5-section massage table.
- The ability to adjust the armrest with one hand - **Smart4Arms**.
- High-quality, soft and pleasant to the touch polyurethane upholstery.
- Progressive foam.
- Two-layer padding (3 cm foam 90 kg/cm<sup>3</sup> + 1 cm foam 43 kg/cm<sup>3</sup>).
- Stitched upholstery.
- Round breathing hole in head section with 2 types of plugs (crescent, round).
- Adjusting the angle of the head section with a gas spring (from -65 to 30°).
- Middle section adjustment angle: -20°.
- Head + middle section adjustment angle: +80°.
- Leg section adjustment angle: +8°.
- Total length of the table 205 cm.
- Maximum safe load of the table when raising-lowering the table ≤ 200 kg.
- Maximum safe working load of the table ≤ 275 kg.
- Safe electric actuator against crushing.
- High/low height adjustment foot bar (49-95 cm).
- F4 wheel system consisting of 4 liftable directional wheels with a central locking mechanism and 4 rubber feet.
- A frame combining wood and steel.
- Wireless remote control with three programmable positions (option).
- Electric adjustment of: head + middle, middle + leg sections via control switches.
- Protection against accidental start-up.
- Linak actuator.

## Settings and dimensions:



## Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Stools for therapist.

## Additional option:

- Wireless remote control.

## Ordering information:



Remote control



Stitched upholstery

**Z5.F4:** 5-section massage table with wheeling system composed of 4 liftable directional wheels with a central locking mechanism and 4 rubber feet.

\* max. safe load when raising-lowering the table



# Theraspace Z7.F4 | Massage table



MD



Box dimensions: **220 x 85 x 80 cm** Weight: **+50 kg**



Length: **205 cm**



Height: **49-95 cm**



Width: **67 cm**



Weight: **110 kg**



Load\*: **200 kg**

## Features:

- 7-section massage table.
- The ability to adjust the armrest with one hand - **Smart4Arms**.
- Swiveling sides in the middle section.
- High-quality, soft and pleasant to the touch polyurethane upholstery.
- Progressive foam.
- Two-layer padding (3 cm foam 90 kg/cm<sup>3</sup> + 1 cm foam 43 kg/cm<sup>3</sup>).
- Stitched upholstery.
- Round breathing hole in head section with 2 types of plugs (crescent, round).
- Adjusting the angle of the head section with a gas spring (from -65 to 30°).
- Middle section adjustment angle: -20°.
- Head + middle section adjustment angle: +80°.
- Leg section adjustment angle: +8°.
- Total length of the table 205 cm.
- Maximum safe load of the table when raising-lowering the table ≤ 200 kg.
- Maximum safe working load of the table ≤ 275 kg.
- Safe electric actuator against crushing.
- High/low height adjustment foot bar (49-95 cm).
- F4 wheel system consisting of 4 liftable directional wheels with a central locking mechanism and 4 rubber feet.
- A frame combining wood and steel.
- Wireless remote control with three programmable positions (option).
- Electric adjustment of: head + middle, middle + leg sections via control switches.
- Protection against accidental start-up.
- Linak actuator.

## Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Stools for therapist.

## Additional option:


- Wireless remote control.

## Ordering information:

**Z7.F4:** 7-section massage table with wheeling system composed of 4 liftable directional wheels with a central locking mechanism and 4 rubber feet.

\* max. safe load when raising-lowering the table

## Settings and dimensions:



Remote control



Stitched upholstery



# Theraspace X4.F4 | Massage table



MD



Box dimensions: **220 x 85 x 80 cm** Weight: **+50 kg**

Length: **205 cm**

Height: **49-95 cm**

Width: **67 cm**


Weight: **110 kg**

Load: **\* 200 kg**

## Features:

- 4-section massage table.
- High-quality, soft and pleasant to the touch polyurethane upholstery.
- Progressive foam.
- Two-layer padding (3 cm foam 90 kg/cm<sup>3</sup> + 1 cm foam 43 kg/cm<sup>3</sup>).
- Stitched upholstery.
- Round breathing hole in head section with 2 types of plugs (crescent, round).
- Adjusting the angle of the head section with a gas spring (from -65 to 30°).
- Head + backrest section adjustment angle: +80°.
- Seat section adjustment angle: +20°.
- Leg section adjustment angle: -10°/+45°.
- Total length of the table 205 cm.
- Maximum safe load of the table when raising-lowering the table ≤ 200 kg.
- Maximum safe working load of the table ≤ 275 kg.
- Safe electric actuator against crushing.
- F4 wheel system consisting of 4 liftable directional wheels with a central locking mechanism and 4 rubber feet.
- High/low height adjustment foot bar (49-95 cm).
- A frame combining wood and steel.
- Wireless remote control with three programmable positions (option).
- Electric adjustment of: head + backrest, seat + leg sections via control switches.
- Protection against accidental start-up.
- Linak actuator.

## Settings and dimensions:



## Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Stools for therapist.

## Additional option:

- Wireless remote control.

## Ordering information:



Remote control



Stitched upholstery

**X4.F4:** 4-section massage table with wheeling system composed of 4 liftable directional wheels with a central locking mechanism and 4 rubber feet.

\* max. safe load when raising-lowering the table



# Theraspace X6.F4 | Massage table



MD



Box dimensions: **220 x 85 x 80 cm** Weight: **+50 kg**



Length: **205 cm**



Height: **49-95 cm**



Width: **67 cm**



Weight: **110 kg**



Load\*: **200 kg**

## Features:

- 6-section massage table.
- The ability to adjust the armrest with one hand - **Smart4Arms**.
- High-quality, soft and pleasant to the touch polyurethane upholstery.
- Progressive foam.
- Two-layer padding (3 cm foam 90 kg/cm<sup>3</sup> + 1 cm foam 43 kg/cm<sup>3</sup>).
- Stitched upholstery.
- Round breathing hole in head section with 2 types of plugs (crescent, round).
- Adjusting the angle of the head section with a gas spring (from -65 to 30°).
- Head + backrest section adjustment angle: +80°.
- Seat section adjustment angle: +20°.
- Leg section adjustment angle: -10°/+45°.
- Total length of the table 205 cm.
- Maximum safe load of the table when raising-lowering the table ≤ 200 kg.
- Maximum safe working load of the table ≤ 275 kg.
- Safe electric actuator against crushing.
- F4 wheel system consisting of 4 liftable directional wheels with a central locking mechanism and 4 rubber feet.
- High/low height adjustment foot bar (49-95 cm).
- A frame combining wood and steel.
- Wireless remote control with three programmable positions (option).
- Electric adjustment of: head + backrest, seat + leg sections via control switches.
- Protection against accidental start-up.
- Linak actuator.

## Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Stools for therapist.

## Additional option:


- Wireless remote control.

## Ordering information:

**X6.F4:** 6-section massage table with wheeling system composed of 4 liftable directional wheels with a central locking mechanism and 4 rubber feet.

\* max. safe load when raising-lowering the table

## Settings and dimensions:



Remote control



Stitched upholstery





# Theraspace X8.F4 | Massage table



MD



Box dimensions: **220 x 85 x 80 cm** Weight: **+50 kg**



Length: **205 cm**



Height: **49-95 cm**



Width: **67 cm**



Weight: **110 kg**




Load\*: **200 kg**

## Features:

- 8-section massage table.
- The ability to adjust the armrests with one hand - **Smart4Arms**.
- Swiveling sides in the backrest section.
- High-quality, soft and pleasant to the touch polyurethane upholstery.
- Progressive foam.
- Two-layer padding (3 cm foam 90 kg/cm<sup>3</sup> + 1 cm foam 43 kg/cm<sup>3</sup>).
- Stitched upholstery.
- Round breathing hole in head section with 2 types of plugs (crescent, round).
- Adjusting the angle of the head section with a gas spring (from -65 to 30°).
- Head + backrest section adjustment angle: +80°.
- Seat section adjustment angle: +20°.
- Leg section adjustment angle: -10°/+45°.
- Total length of the table 205 cm.
- Maximum safe load of the table when raising-lowering the table ≤ 200 kg.
- Maximum safe working load of the table ≤ 275 kg.
- Safe electric actuator against crushing.
- F4 wheel system consisting of 4 liftable directional wheels with a central locking mechanism and 4 rubber feet.
- High/low height adjustment foot bar (49-95 cm).
- A frame combining wood and steel.
- Wireless remote control with three programmable positions (option).
- Electric adjustment of: head + backrest, seat + leg sections via control switches.
- Protection against accidental start-up.
- Linak actuator.

## Settings and dimensions:



## Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Stools for therapist.

## Additional option:

- Wireless remote control.

## Ordering information:



Remote control



Stitched upholstery

**X8.F4:** 8-section massage table with wheeling system composed of 4 liftable directional wheels with a central locking mechanism and 4 rubber feet.


\* max. safe load when raising-lowering the table





# Theraspace

## Line



### Technical data:

	E-S2	M-S2	M-S4	Z3	Z5	Z7	X4	X6	X8
Dimensions (length x width) [cm]:	202 x 69	203 x 69	203 x 69	205 x 67	205 x 67	205 x 67	205 x 67	205 x 67	205 x 67
Number of sections:	2	2	4	3	5	7	4	6	8
Head section adjustment angle [°]:	-20/+85	-70/+40	-70/+40	-65/+30	-65/+30	-65/+30	-65/+30	-65/+30	-65/+30
Middle section adjustment angle [°]:	-	-	-	-20	-20	-20	-	-	-
Head + middle s. adjustment angle [°]:	-	-	-	+80	+80	+80	-	-	-
Head + backrest s. adjustment angle [°]:	-	-	-	-	-	-	+80	+80	+80
Seat section adjustment angle [°]:	-	-	-	-	-	-	+20	+20	+20
Leg section adjustment angle [°]:	-	-	-	+8	+8	+8	-10/+45	-10/+45	-10/+45
Height adjustment [cm]:	47-95	47-95	47-95	49-95	49-95	49-95	49-95	49-95	49-95
Weight [kg]:	115	115	115	110	110	110	110	110	110
Max. lifting capacity [kg]:	200	200	200	200	200	200	200	200	200
Power supply [V/Hz]:	100-240/50-60	100-240/50-60	100-240/50-60	100-240/50-60	100-240/50-60	100-240/50-60	100-240/50-60	100-240/50-60	100-240/50-60
Max. working load [kg]:	275	275	275	275	275	275	275	275	275

### Upholstery OMEGA:

Surface composition: 100% Polyurethane      Backing composition: 100% Polyester

OM006 white	OM240 beige	OM315 sky	OM221 blue	OM314 cement	OM076 olive grey	OM773 dark brown	OM225 black	OM469 lima	OM157 forest	OM053 red	OM217 bordeaux
----------------	----------------	--------------	---------------	-----------------	---------------------	---------------------	----------------	---------------	-----------------	--------------	-------------------



**Durability and performance:**  
>500 000 cycles (intensive usage)



**Resistance to hydrolysis** (8 years)



**High scratch resistance**



**Does not contain chemical solvents**




**Respects the environment**

## Types of wheeling system:


### F0

F0 system composed of 2 legs with rubber feet and 2 small wheels for moving of Safari Table by lifting it from one side. The system is recommended for users who relocate their Safari Table less frequently.




## Massage tables

### Safari




## F4

System composed of 4 directional castors with central brake allows you very easy and fast relocation of Safari Table with or without a patient.



2 side brackets for fixation of stabilization belts and elastic bands

Two-layer padding with upholstery in many colours



Foot - operated lever of F4 system


Rubber foot with 1 cm height adjustment



# Safari Puma S2.F0/S2.F4

Massage table

MD



Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**



Length: **202 cm**



Height: **52-93 cm**



Width: **69 cm**



Weight: **86/99 kg**



Max load\*: **200 kg**

## Features:

- 2-section massage table.
- Width of table top: 69 cm.
- Adjustable head section from -70° to +40° with a gas spring.
- Round face hole in the head section with a removable face hole cover which increases patient's comfort (2 types of cover: round and crescent).
- Two-layer padding with upholstery in many colours.
- Powder coated, stable scissor steel frame.
- 2 side brackets for the fixation of stabilization belts and elastic bands
- Traction support sleeves.
- Personal Authorization Safety System with 2 magnetic safety keys.
- HalloTronic® System electric high/low height adjustment from 52 to 93 cm with a lifting capacity up to 200 kg.
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table.
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.
- Two ways of frame painting: „white“ or „twintone“.


## Additional option:

- Built-in surface heating system of table top.
- Upholstery of 80 mm thickness.
- Upholstery of 80 cm width.

## Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Traction supports.
- Stabilization belts.
- Stools for therapist.

## Settings and dimensions:



## Ordering information:

**S2.F0:** 2-section massage table with wheeling system composed of 2 small wheels and 2 legs with rubber feet.

**S2.F4:** 2-section massage table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.

\* max. lifting capacity



# Safari

## Puma S4.F0/S4.F4

Massage table

MD



Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**



Length: **202 cm**



Height: **52-93 cm**



Width: **69 cm**



Weight: **89/102 kg**



Max load\*: **200 kg**

### Features:

- 4-section massage table.
- Width of table top: 69 cm.
- Adjustable head section from -70° to +40° with a gas spring.
- Round face hole in the head section with a removable face hole cover which increases patient's comfort (2 types of cover: round and crescent).
- Rotating slide-down armrest in head section.
- Two-layer padding with upholstery in many colours.
- Powder coated, very stable scissor steel frame.
- 2 side brackets for the fixation of stabilization belts and elastic bands.
- Traction support sleeves.
- Personal Authorization Safety System with 2 magnetic safety keys.
- HalloTronic® System electric high/low height adjustment from 52 to 93 cm with a lifting capacity up to 200 kg.
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table.
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.
- Two ways of frame painting: „white“ or „twintone“.

### Additional option:

- Built-in surface heating system of table top.
- Upholstery of 80 mm thickness.
- Upholstery of 80 cm width.

### Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Traction supports.
- Stabilization belts.
- Stools for therapist.


### Ordering information:

**S4.F0:** 4-section massage table with wheeling system composed of 2 small wheels and 2 legs with rubber feet.

**S4.F4:** 4-section massage table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.

\* max. lifting capacity

### Settings and dimensions:





# Safari

## Puma E-S2.F0/E-S2.F4

Massage table

MD



2 ways of frame painting



Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**



Length: **208 cm**



Height: **52-93 cm**



Width: **69 cm**



Weight: **88/103 kg**



Max load\*: **200 kg**

### Features:

- 2-section massage table with a long head section.
- Width of table top: 69 cm.
- Adjustable long head section from -20° to +85° or 0° to +85° (with special frame) with 2 gas springs.
- Round face hole in the head section with a removable face hole cover (2 types: round and crescent).
- Two-layer padding with upholstery in many colours.
- Powder coated, extremely stable steel frame.
- 2 side brackets for the fixation of stabilization belts and elastic bands.
- Personal Authorization Safety System with 2 safety magnetic keys.
- HalloTronic® System for electric high/low height adjustment from 52 to 93 cm with a lifting capacity up to 200 kg.
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table.
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.
- Two ways of frame painting: „white" or „twintone".


### Additional option:

- Upholstery of 80 mm thickness.
- Upholstery of 80 cm width.
- Built-in adjustable heating system.

### Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Traction supports.
- Stabilization belts.
- Stools for therapist.

### Settings and dimensions:



### Ordering information:

**E-S2.F0:** 2-section massage table with wheeling system composed of 2 small wheels and 2 legs with rubber feet.

**E-S2.F4:** 2-section massage table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.

\* max. lifting capacity





# Safari

## Leopard S3.F0/S3.F4

Massage table



MD



2 ways of frame painting



Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**



Length: **205 cm**



Height: **52-93 cm**



Width: **69 cm**



Weight: **93/106 kg**



Max load\*: **200 kg**

### Features:

- 3-section massage table.
- Width of table top: 69 cm.
- Adjustable head section from -70° to +40° with a gas spring.
- Round face hole in the head section with a removable face hole cover which increases patient's comfort (2 types of cover: round and crescent).
- Adjustable leg section from 0° to 90° with a gas spring.
- Two-layer padding with upholstery in many colours.
- Powder coated, very stable scissor steel frame.
- 2 side brackets for the fixation of stabilization belts and elastic bands.
- Traction support sleeves.
- Personal Authorization Safety System with 2 magnetic safety keys.
- HalloTronic® System electric high/low height adjustment from 52 to 93 cm with a lifting capacity up to 200 kg.
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table.
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.
- Two ways of frame painting: „white" or „twintone".

### Additional option:

- Built-in surface heating system of table top.
- Upholstery of 80 mm thickness.
- Upholstery of 80 cm width.

### Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Traction supports.
- Stabilization belts.
- Stools for therapist.


### Ordering information:

**S3.F0:** 3-section massage table with wheeling system composed of 2 small wheels and 2 legs with rubber feet.

**S3.F4:** 3-section massage table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.

\* max. lifting capacity

### Settings and dimensions:



# Safari

## Leopard S5.F0/S5.F4

Massage table



MD



2 ways of frame painting



Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**



Length: **205 cm**



Height: **52-93 cm**



Width: **69 cm**



Weight: **96/109 kg**



Max load\*: **200 kg**

### Features:

- 5-section massage table.
- Width of table top: 69 cm.
- Adjustable head section from -70° to +40° with a gas spring.
- Round face hole in the head section with a removable face hole cover which increases patient's comfort (2 types of cover: round and crescent).
- Rotating slide-down armrests.
- Adjustable leg section from 0° to 90° with a gas spring.
- Two-layer padding with upholstery in many colours.
- Powder coated, very stable scissor steel frame.
- 2 side brackets for the fixation of stabilization belts and elastic bands.
- Traction support sleeves.
- Personal Authorization Safety System with 2 magnetic safety keys.
- HalloTronic® System electric high/low height adjustment from 52 to 93 cm with a lifting capacity up to 200 kg.
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table.
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.
- Two ways of frame painting: „white“ or „twintone“.

### Additional option:

- Built-in surface heating system of table top.
- Upholstery of 80 mm thickness.
- Upholstery of 80 cm width.

### Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Traction supports.
- Stabilization belts.
- Stools for therapist.


### Ordering information:

**S5.F0:** 5-section massage table with wheeling system composed of 2 small wheels and 2 legs with rubber feet.

**S5.F4:** 5-section massage table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.

\* max. lifting capacity

### Settings and dimensions:



# Safari

## Jaguar P3.F0/P3.F4

Massage table



MD



2 ways of frame painting



Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**



Length: **205 cm**



Height: **54-95 cm**



Width: **69 cm**



Weight: **99/112 kg**



Max load\*: **200 kg**

### Features:

- 3-section massage table with postural drainage setting.
- Postural drainage with an electric actuator operated by 2 switches from both sides of the table.
- Width of table top: 69 cm.
- Adjustable head section from -70° to +40° with a gas spring.
- Round face hole in the head section with a removable face hole cover which increases patient's comfort (2 types of cover: round and crescent).
- Adjustable leg section from 0° to 90° with gas spring.
- Two-layer padding with upholstery in many colours.
- Powder coated, very stable scissor steel frame.
- 2 side brackets for the fixation of stabilization belts and elastic bands.
- Traction support sleeves.
- Personal Authorization Safety System with 2 magnetic safety keys.
- HalloTronic® System electric high/low height adjustment from 54 to 95 cm with a lifting capacity up to 200 kg.
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table.
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.
- Two ways of frame painting: „white“ or „twintone“.

### Additional option:

- Built-in surface heating system of table top.
- Upholstery of 80 mm thickness.
- Upholstery of 80 cm width.

### Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Traction supports.
- Stabilization belts.
- Stools for therapist.


### Ordering information:

**P3.F0:** 3-section massage table with wheeling system composed of 2 small wheels and 2 legs with rubber feet.

**P3.F4:** 3-section massage table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.

\* max. lifting capacity

### Table articulations:



# Safari

## Jaguar P5.F0/P5.F4

Massage table

MD



2 ways of frame painting



Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**



Length: **205 cm**



Height: **54-95 cm**



Width: **69 cm**



Weight: **102/115 kg**



Max load\*: **200 kg**

### Features:

- 5-section massage table with the postural drainage setting.
- Postural drainage with an electric actuator operated by 2 switches from both sides of the table.
- Width of table top: 69 cm.
- Adjustable head section from -70° to +40° with a gas spring.
- Round face hole in the head section with a removable face hole cover which increases patient's comfort (2 types of cover: round and crescent).
- Rotating slide-down armrests.
- Adjustable leg section from 0° to 90° with a gas spring.
- Two-layer padding with upholstery in many colours.
- Powder coated, very stable scissor steel frame.
- 2 side brackets for the fixation of stabilization belts and elastic bands.
- Traction support sleeves.
- Personal Authorization Safety System with 2 safety magnetic keys.
- HalloTronic® System electric high/low height adjustment from 54 to 95 cm with a lifting capacity up to 200 kg.
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table.
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.
- Two ways of frame painting: „white“ or „twintone“.

### Additional option:

- Built-in surface heating system of table top.
- Upholstery of 80 mm thickness.
- Upholstery of 80 cm width.

### Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Traction supports.
- Stabilization belts.
- Stools for therapist.


### Ordering information:

**P5.F0:** 5-section massage table with wheeling system composed of 2 small wheels and 2 legs with rubber feet.

**P5.F4:** 5-section massage table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.

\* max. lifting capacity

### Table articulations:



# Safari

## Line




### Technical data:

	Puma S2	Puma E-S2	Puma S4	Leopard S3	Leopard S5	Jaguar P3	Jaguar P5
Dimensions (length x width) [cm]:	202 x 69	208 x 69	202 x 69	205 x 69	205 x 69	205 x 69	205 x 69
Number of sections:	2	2	4	3	5	3	5
Head section inclination angle [°]:	-70/+40	-20/+85	-70/+40	-70/+40	-70/+40	-70/+40	-70/+40
Middle section inclination angle [°]:	-	-	-	-	--	0/+24	0/+24
Leg section inclination angle [°]:	-	-	-	0/+90	0/+90	0/+90	0/+90
Height adjustment [cm]:	52-93	52-93	52-93	52-93	52-93	54-95	54-95
Weight F0/F4 [kg]:	86/99	88/103	89/102	93/106	96/109	99/112	102/115
Max. lifting capacity [kg]:	200	200	200	200	200	200	200
Power supply [V/Hz]:	100-240/50-60	100-240/50-60	100-240/50-60	100-240/50-60	100-240/50-60	100-240/50-60	100-240/50-60
Max. working load [kg]:	275	275	275	275	275	275	275



Adjustment of head rest



Adjustment of postural drainage



2 types of breathing slot cover:  
round and crescent



Built-in surface heating system  
of table top (option)



Personal Authorization  
Safety System (PASS)



# Safari Puma/Leopard/Jaguar

Massage table  
Hydraulic



MD



2 ways of frame painting



Box dimensions: **220 x 85 x 80 cm** Weight: **+50 kg**



Length: **202-208 cm**



Height: **52-99 cm / 54-101 cm**



Width: **69 cm**



Weight: **92-121 kg**




Max load\*: **200 kg**


## Features:

- 2-, 3-, 4- or 5-section massage table with a powerful hydraulic actuator.
- Postural drainage by the means of 2 levers accessible from 2 sides of table top (P3, P5 version).
- Rotating armrests in head section (S4, S5, P5 versions).
- Width of table top: 69 cm.
- Adjustable head section from -70° to +40° with a gas spring (-20° to 85° in E-S2 version (tables with long head section)).
- Round face hole in the head section with a removable face hole cover which increases patient's comfort (2 types of cover: round and crescent).
- Adjustable leg section from 0° to 90° with a gas spring.
- Two-layer padding with upholstery in many colours.
- Powder coated, very stable scissor steel frame.
- 2 side brackets for the fixation of stabilization belts and elastic bands.
- Traction support sleeves.
- Hydraulic high/low height adjustment with a lifting capacity up to 200 kg.
- High/low height adjustment with foot pedal accessible from both sides of the table.
- Two ways of frame painting: „white” or „tintone”.

## Frame F0-H:



## Frame F4-H:



## Additional option:

- Upholstery of 80 mm thickness.
- Upholstery of 80 cm width.

## Technical data:

	Puma S2	Puma E-S2	Puma S4	Leopard S3	Leopard S5	Jaguar P3	Jaguar P5
Dimensions (length x width) [cm]:	202 x 69	208 x 69	202 x 69	205 x 69	205 x 69	205 x 69	205 x 69
Number of sections:	2	2	4	3	5	3	5
Head section inclination angle [°]:	-70/+40	-20/+85	-70/+40	-70/+40	-70/+40	-70/+40	-70/+40
Middle section inclination angle [°]:	-	-	-	-	-	0/+24	0/+24
Leg section inclination angle [°]:	-	-	-	0/+90	0/+90	0/+90	0/+90
Height adjustment [cm]:	52-99	52-99	52-99	52-99	52-99	54-101	54-101
Weight F0/F4 [kg]:	92/105	96/109	95/108	99/112	102/115	105/118	108/121
Max. lifting capacity [kg]:	200	200	200	200	200	200	200

## Ordering information:

**S2.F0-H**  
**S2.F4-H**

**E-S2.F0-H**  
**E-S2.F4-H**

**S4.F0-H**  
**S4.F4-H**

**S3.F0-H**  
**S3.F4-H**

**S5.F0-H**  
**S5.F4-H**

**P3.F0-H**  
**P3.F4-H**

**P5.F0-H**  
**P5.F4-H**

\* max. lifting capacity





# WaveMotion Ellipse

Massage table



MD





Box dimensions: **220 x 95 x 95 cm** Weight: **+50 kg**

Massage table **WaveMotion Ellipse** - gentle swing and rotation 360° in space. This massage table is designed for performing manual massage treatments to treat and alleviate traumas or disabilities by restoring or improving the function of massaged tissues, organs, and systems. It enables the patient to take the correct position for massage treatment and to adjust height and/or position of the table sections.

Manual massage treatments using this table has a positive effect on i.a.:

- muscular tension disorders,
- spinal disc disorder,
- chronic back, low back pain,
- cystic fibrosis,
- pain in the quadriceps fatigue,
- tension-type headache,
- osteoarthritis of the hip joints.



## Features:

- Versions depended on the therapist needs: 2 - or 4 - sections, upholstery padding of 4 or 8 cm thickness, width 75 cm or 80 cm (option).
- Easily adjustable sections: head section - via gas spring, arm rests – via rotating mechanism.
- Table top movement (360° rotation on a floating suspension and can be tilt both axes).
- Patient comfort: round face hole with a matching cushion.
- Two-layer padding (specially selected density of foams), with upholstery in many colours.
- Stable steel frame.
- Personal Authorization Safety System (option).
- Electric height adjustment via foot controller.

## Technical data:

	S2	S4
Dimensions (L x W) [cm]:	202 x 75	202 x 75
Section number:	2	4
Fold-down armrests:	-	+
Head section adjustment:	+	+
Height adjustment [cm]:	65-105	65-105
Max. weight [kg]:	145	145
Max. lifting capacity [kg]:	150	150
Power supply [V/Hz]:	230/50/60	230/50/60

## Models available:

- **S2:** 2-section with adjustable head section,
- **S4:** 4-section with adjustable head section and fold-down armrests.

**S2, S4 models:** the 40°/90° headrest allows for an ergonomic positioning in the prone position, but the head of patient/client can be gently supported in the supine position. The adjustment range of the headrest is from 40° positive to 90° negative. The headrest has a round face cradle with a matching cushion.

## Additional options:

- Upholstery resistant to oils (Premium Upholstery).
- Built-in surface heating system of table top.
- Upholstery of 80 cm width.
- F2 lift (2 fixed wheels + 2 castors central lifted).
- Personal Authorization Safety System.
- Additional foot controller.



# Nexus


## S2/E-S2/S4/S2-H/E-S2-H/S4-H

Massage table

MD



Electric height adjustment



Hydraulic height adjustment

Massage table **Nexus** is designed for performing manual, massage treatments to treat and alleviate traumas or disabilities by restoring or improving the function of massaged tissues.

### Features:

- 2- or 4-section massage table.
- Available in electric or hydraulic version.
- Electric height adjustment with a hand - or foot-operated controller or hydraulic height adjustment with a lever.
- Safe electric Linak actuator against crushing (only in electric version).
- Head section adjustment with a gas spring.
- Rotating slide-down armrests in 4-section table.
- Stable scissor base.
- Durable construction with lifting capacity up to 185 kg.
- Mobile - 2 wheels.
- Wide range of upholstery colours.

### Technical data:

	S2/S2-H	E-S2/E-S2-H	S4/S4-H
Number of sections:	2	2	4
Head section angle adjustment [°]:	-70/+40	-10/+85	-70/+40
Height adjustment (el./H) [cm]:	50-94	50-94	50-94
Dimensions (L x W) [cm]:	200 x 69	210 x 69	200 x 69
Weight [kg]:	85	85	85
Max. lifting capacity [kg]:	185	185	185
Power supply [V/Hz] (only electric version):	100-240/50-60	100-240/50-60	100-240/50-60


### Additional option:

- Built-in surface heating system of table top (available only in electric version).


### Ordering information:

### Table articulation:

#### Nexus type S2/S2-H



#### Nexus type S4/S4-H



#### Nexus type E-S2/E-S2-H



### Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Stools for therapist.
- Optional: 2 side brackets or 2 horizontal folding barriers.



Horizontal folding barrier

#### S2:

2-section massage table, electric height adjustment.

#### S2-H:

2-section massage table, hydraulic height adjustment.

#### E-S2:

2-section massage table (long head section), electric height adjustment.

#### E-S2-H:

2-section massage table (long head section), hydraulic height adjustment.

#### S4:

4-section massage table, electric height adjustment.

#### S4-H:

4-section massage table, hydraulic height adjustment.

\* max. lifting capacity



# Safari Elephant B-S1.F0/ B-S1.F4/B-S1.F0-H/B-S1.F4-H

Neurological  
rehabilitation table

MD



Electric height adjustment



Hydraulic height adjustment



Box dimensions: **220 x 136 x 75 cm** Weight: **+90 kg**



Length: **200 cm**



Height: **50-99 cm**



Width: **120 cm**



Weight: **150 kg**



Max load\*: **200 kg**

Neurological rehabilitation table **Safari type Elephant** is designed for improving the motor skills of children with central nervous system impairment as well as Bobath and Vojta treatment.

## Features:

- 1-section table (electric or hydraulic).
- Extra wide surface.
- Two-layer padding with upholstery in many colours.
- Powder coated, extremely wide and stable steel frame.
- Personal Authorization Safety System with 2 magnetic safety keys (only in el. version).
- HalloTronic® System for electric high/low height adjustment from 50 to 99 cm (el. version).
- Hydraulic high/low height adjustment (50-99 cm) via foot pedal accessible from both sides of the table (hydr. version).
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table (el. version).
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.
- F4 wheeling system - dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.
- Max. lifting capacity up to 200 kg.
- Two ways of frame painting: „white“ or „twintone“.
- Max. working load up to 275 kg.

## Additional option:

- Face hole.
- Built-in surface heating system of table top (only in el. version).
- Upholstery of 80 mm thickness.


## Ordering information:

**B-S1.F0/B-S1.F0-H:** 1-section electric/hydraulic neurological rehabilitation table with basic wheeling system composed of 2 small wheels and 2 legs with rubber feet.

**B-S1.F4/ B-S1.F4-H:** 1-section electric/hydraulic neurological rehabilitation table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.

\* max. lifting capacity

## Table articulations:



## Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Stools for therapist.



# Safari Elephant B-S2.F0/ B-S2.F4/B-S2.F0-H/B-S2.F4-H

Neurological  
rehabilitation table

MD



Electric height adjustment



Hydraulic height adjustment



Box dimensions: **220 x 136 x 75 cm** Weight: **+90 kg**



Length: **200 cm**



Height: **50-99 cm**



Width: **120 cm**



Weight: **150 kg**



Max load\*: **200 kg**

Neurological rehabilitation table **Safari type Elephant** is designed for improving the motor skills of children with central nervous system impairment as well as Bobath and Vojta treatment.

## Features:

- 2-section table (electric or hydraulic).
- Extra wide surface.
- Adjustable head section from 0° to +85° via 2 gas springs.
- Two-layer padding with non-flammable, bio-compatible and scratchproof upholstery in many colours.
- Powder coated, extremely wide and stable steel frame.
- Personal Authorization Safety System with 2 magnetic safety keys (only in el. version).
- HalloTronic® System for electric high/low height adjustment from 50 to 99 cm (electric version).
- Hydraulic high/low height adjustment (50-99 cm) via foot pedal accessible from both sides of the table (hydraulic version).
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table (electric version).
- Rubber feet with 1 cm height adjustment allowing for the levelling of the table.
- F4 wheeling system - dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.
- Max. lifting capacity up to 200 kg.
- Two ways of frame painting: „white” or „twintone”.


## Additional option:

- Face hole.
- Built-in surface heating system of table top (only in el. version).
- Upholstery of 80 mm thickness.

## Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Stools for therapist.

## Table articulations:







## Ordering information:

- B-S2.F0/B-S2.F0-H:** 2-section electric/hydraulic neurological rehabilitation table with basic wheeling system composed of 2 small wheels and 2 legs with rubber feet.
- B-S2.F4/B-S2.F4-H:** 2-section electric/hydraulic neurological rehabilitation table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with rubber feet.

\* max. lifting capacity



MD



7-section chiropractic table **Chirospace** is a multifunctional equipment that ensures comfort of work for every therapist. Thanks to the ergonomic construction of the table top, with the possibility of lateral and vertical rotation of the foot section, the physiotherapist can adjust the initial correction of the patient's spine, work on soft tissues, and perform any manual procedures.

## Features:

- 7-section table.
- Head section divided into 2 parts (adjustable gap: 2-6 cm).
- 2 possibilities of head section adjustment: from -35° to +25° with a gas spring and additional adjustment from 0° to +18° (in opposite direction) via second gas spring.
- Rotating slide down armrests (adjustment range: up to -24 cm vertical).
- 2 possibilities of leg section adjustment: from 0° to -25° and side rotation from -9° to 9°.
- Electric high/low height adjustment via foot-operated controller at two sides of the table base.
- 3-layer padding 40 mm, with upholstery in many colours.
- Shift of shank support: 0-19 cm.
- Settings lock: rotation and shift.
- DROP option (2 cm).
- Second DROP option.
- Mechanism of dynamic force compensation - to minimize the patient's body weight resting on the leg section of the table via an offload and steering mechanism.
- 2 fixed wheels + 2 castors central lifted (F2 system).

## Technical data:

Dimensions (LxW) [cm]:	185-210 x 50 (55 - head section + armrest)
Height adjustment [cm]:	55-96
Weight [kg]:	135
Max. lifting capacity [kg]:	180
Shift [cm]:	0-19
Power supply [V; Hz]:	100-240; 50-60

\* max. lifting capacity



## Additional accessories:

- Wedges.
- Bolsters, semi-bolsters.
- Stools for therapist.



# Vertimo Classic/Hi-Lo | Tilting tables



MD



1-section (UNO)

2-section (DUO)

2-section (DUO+)



Box dimensions: **220 x 85 x 80 cm** Weight: **+50 kg**

Tilting table **Vertimo type Classic** - very stable 1-section tilting table with electric adjustment of verticalisation.

## Features:


- 1-section table.
- Electric actuator for tilting adjustment.
- Manual remote control for the adjustment of tilting angle with access authorization - with a special key placed in the remote control.
- Brackets for the attachment of fastening belts (standard accessories).
- Two-layer padding with upholstery in many colours.
- Powder coated, very stable steel frame.
- 4 large separately lockable castors (rolling and rotation lock).
- Emergency power supply battery.

## Technical data:

Footrest inclination angle [°]:	-30/+15
Number of sections:	1
Dimensions (L x W) [cm]:	190 x 69
Table top inclination angle [°]:	0-87 (+3)
Height [cm]:	53 (fixed)
Max. table load [kg]:	150
Weight [kg]:	100
Power supply [V/Hz]:	230/50/60

## Additional accessories for Vertimo line:

- Adjustable support under the arm.
- Top of the small table.
- Stools for therapist.



Small table



Arm supports



Box dimensions: **220 x 85 x 80 cm** Weight: **+50 kg**

Tilting table **Vertimo type HI-LO** - very stable 1- or 2-section tilting table with electric adjustment of verticalization. Available in 3 versions:

- 1-section (**UNO**) tilting table,
- 2-section (**DUO**) tilting table,
- 2-section (**DUO+**) tilting table.

## Features:

- 1 or 2-section table.
- 2 electric actuators – one ensures smooth tilting; the other one adjusts height.
- Electric height and tilting adjustment with manual remote control with the function of access authorization - authorization is done by a special key placed in the remote control.
- Two-layer padding with upholstery in many colours.
- Powder coated, very stable steel frame.
- 2 separately adjustable footrests.
- Brackets for the attachment of fastening belts (standard accessories).
- 4 large separately lockable castors (rolling and rotation lock).
- Tilt angle indicator.
- Emergency power supply battery (allows table lowering to min. height and return to horizontal position at power failures).
- Vertimo Hi-Lo DUO+ - table with faster actuator, quick release system (ensures fast table top lowering in emergency situations), bracket for leg belts.

## Technical data:

Number of sections:	1 or 2
Dimensions (L x W) [cm]:	190 x 69
Footrest length [cm]:	36
Table top inclination angle [°]:	0-87 (+3)
Head section inclination angle (2-section ver.) [°]:	0-85
Footrest inclination angle [°]:	-30/+15
Height adjustment (1-, 2-section ver./DUO+ ver.) [cm]:	49,5-96,5/54-92
Max. table load [kg]:	150
Castors diameter [cm]:	10
Weight [kg]:	130
Power supply [V/Hz]:	230/50/60



# Fiji P5/P7

Pacific examination  
and treatment tables

# Hawaii S2/S4

Pacific examination  
and treatment tables



MD



Box dimensions: **220 x 95 x 95 cm** Weight: **+50 kg**



Length: **212 cm**



Height: **64-104 cm**



Width: **80 cm**



Weight: **150/130 kg**



Max load\*: **150 kg**


**Pacific** examination and treatment table is intended for standard examinations and procedures carried out in a wide scope, also used for simple medical activities. It is usually used in professional medical care facilities, i.e. hospital, clinic, as well as in doctor's offices, treatment rooms or diagnostic laboratories.

## Additional option:

- Integrated option of heating and LED lights.
- Central lift (2 castors central lifted).

## Features:

- 5- or 7-section examination and treatment table with electric height adjustment by the means of a remote (controller).
- Head section adjustment by the means of gas spring.
- Head rest with shaped face hole with crescent plug.
- Foldable armrests lowered till 19 cm (only in 7 section table) .
- Adjustment of: back, thighs, shanks sections by means of electric actuators.
- Multi-layered and very soft padding of 15 cm thickness with upholstery in many colours.
- Trendelenburg/reverse Trendelenburg position (15°).
- 2 wheels.




P7

## Technical data:

	FIJI	HAWAII
Number of sections:	5 or 7	2 or 4
Height adjustment [cm]:	64-104	64-104
Dimensions (L x W) [cm]:	212 x 80	212 x 80
Weight [kg]:	150	130
Max. lifting capacity [kg]:	150	150
Power supply [V/Hz]:	120-240 / 50/60	120-240 / 50/60

## Features:

- 2-section or 4-section examination and treatment table with electric height adjustment by the means of a remote (controller).
- Head section adjustment by the means of gas spring.
- Head rest with shaped face hole with crescent plug.
- Foldable armrests lowered till max. 19 cm (only in version 4-section).
- Multi-layered and very soft padding of 15 cm thickness with upholstery in many colours.
- Trendelenburg/reverse Trendelenburg position (15°).
- 2 wheels.




P4

\* max. lifting capacity



# Kangoo S3.C4 | Examination and treatment table

MD



Box dimensions: **220 x 85 x 75 cm** Weight: **+50 kg**



Length: **189 cm**



Height: **46-91 cm**



Width: **72 cm**



Weight: **116 kg**



Max load\*: **250 kg**

Examination and treatment table **Kangoo model S3.C4** is intended for standard examinations and procedures carried out in a wide scope. Well-thought-out design, the possibility to adjust the angles of individual table sections and the stability of the product ensure that the patient is in a comfortable position for examination/treatment.

The table offers possibility of the Trendelenburg and reverse Trendelenburg position (18°).


## Features:

- 3- section table, with a long head section.
- Adjustable head and foot section from -18° to 85° with a gas springs.
- Trendelenburg/reverse Trendelenburg.
- Face hole in the head section with a removable face hole cover.
- Two-layer padding with upholstery in many colours.
- Powder coated stable steel frame.
- 2 folding horizontal barriers on both sides; fold-down when table subjected to normal use activities.
- Electric height adjustment (47-92 cm) using manual remote control.
- Battery backup, cordless operation.
- Safety - table equipped with an access lock - safety key and automatic locking ( to lock height adjustment function).
- 4 large 125 mm braked castors with foot locking mechanism.

## Technical data:

Height adjustment [cm]:	46-91
Dimensions (length x width) [cm]:	189 x 72
Weight [kg]:	116
Max. lifting capacity [kg]:	250
Paper roll width [cm]:	60
Power supply [V/Hz]:	100-240/50-60

\* max. lifting capacity



## Additional accessories:

- Paper holder.
- Wedges.
- Bolsters, semi-bolsters.
- Stools for therapist.
- Basket e.g. for clothes.



# Safari Puma/Leopard

Treatment and examination table

MD



Box dimensions: **220 x 85 x 80 cm** Weight: **+50 kg**



Length: **204-208 cm**



Height: **52-93 cm**



Width: **69 cm**



Weight: **≤ 115 kg**



Max load\*: **200 kg**

## Features:


- 2-, 3-, 4- or 5-section table.
- Rotating armrests in head section (S4, S5 versions).
- Width of table top: 69 cm.
- Adjustable head section from -70° to +40° with a gas spring (-20° to 85° in E-S2 version (tables with long head section)).
- Round face hole in the head section with a removable face hole cover which increases patient's comfort (2 types of cover: round and crescent).
- Adjustable leg section from 0° to 90° with a gas spring.
- Two-layer padding with upholstery in many colours.
- Powder coated, very stable scissor steel frame.
- 2 side brackets for the fixation of stabilization belts and elastic bands.
- Traction support sleeves.
- Electric high/low height adjustment with a lifting capacity up to 200 kg.
- At bumpers option head section can't be lower below 0°.

## Technical data:


	Puma S2	Puma E-S2	Puma S4	Leopard S3	Leopard S5
Dimensions (length x width) [cm]:	204 x 69	208 x 69	204 x 69	205 x 69	205 x 69
Number of sections:	2	2	4	3	5
Head section inclination angle [°]:	-70/+40	-20/+85	-70/+40	-70/+40	-70/+40
Middle section inclination angle [°]:	-	-	-	-	-
Leg section inclination angle [°]:	-	-	-	0/+90	0/+90
Height adjustment [cm]:	52-93	52-93	52-93	52-93	52-93
Weight [kg]:	≤ 115	≤ 115	≤ 115	≤ 115	≤ 115
Max. lifting capacity [kg]:	200	200	200	200	200
Power supply [V/Hz/A]:	100-240/50-60/3-1,5				

## Additional option:

- Upholstery of 80 mm thickness.
- Upholstery of 80 cm width.
- Head section without face hole.
- Wheels of bigger diameter (C version).
- Hand- or foot- operated controller.
- Bumpers.
- Equipotential pin.
- CPR (for lowering the table enabling proper cardiopulmonary resuscitation).



Equipotential pin



CPR

## Additional accessories:

- Bolsters, wedges.
- Traction supports.
- Horizontal folding barriers.
- Rail for mounting accessories.
- Paper holder.
- Stabilization belts.

## Ordering information:

**S2.M4**   **E-S2.M4**   **S4.M4**   **S3.M4**   **S5.M4**  
**S2.M4-C**   **E-S2.M4-C**   **S4.M4-C**   **S3.M4-C**   **S5.M4-C**

\* max. lifting capacity







# Safari Buffalo E2.C4-H

Patient transport table



MD



Patient transport table **Safari type Buffalo** is designed for efficient transport in medical facilities. The table ensures easy operation for caregivers. The simple design provides reduced effort for medical personnel, improved quality of patient care and greater comfort for the patient.

## Features:

- 2-sections transport table with a powerful hydraulic actuator.
- Adjustable, long head section from 0° to +70° with 2 gas springs .
- Removable mattress (60 mm thickness) made of high-quality polyurethane foam (40-45 kg/m³) and fabric in black colour; fabric features: flame retardant/mildew resistant/impermeable/easy to clean (can be washed in 60-90°C, resistant for disinfectants).
- Powder coated, very stable scissor steel frame.
- Hydraulic high/low height adjustment with a lifting capacity up to 200 kg.
- High/low height adjustment from 51 to 100 cm with foot pedal accessible from both sides of the table.
- 4 large 125 mm TENTE noiseless castors (total lock) which increase manoeuvrability; 1 electro-conductive (anti-static) wheel marked with a yellow-coloured border, to prevent electrostatic discharge.

## Accessories:



Hanger for oxygen tank



Basket for i.e. clothes




Drip hanger (can be attached on both sides in easy and fast way)



600 mm wide paper towel hanger located on the head section

\* max. lifting capacity

## Table articulations:



## Equipped with:

- Folding horizontal barriers on both sides
- Handles for transport on the head and leg sections
- Bumpers at the corners

## Additional option:

- The table offers possibility of the Trendelenburg/reverse Trendelenburg position (-10°/+10°).

## Technical data:

Dimensions (length x width) [cm]:	228 x 80
Number of sections:	2
Head section inclination angle [°]:	0/+70
Laying width [cm]:	69
Height adjustment [cm]:	51-100
Weight [kg]:	135
Max. lifting capacity [kg]:	200



MD



Box dimen.: **198 x 66 x 17 cm** Weight: **+2 kg**

## LZM type 3:

**LZM** couches are designed to be used in therapeutic procedures to heal or provide relief in injuries or other impairments by means of natural physical factors. The couch is also used in doctor's surgeries. It allows the patients to assume a comfortable horizontal position thanks to its adjustable head rest.

## Features:

- Adjustable head rest (10 steps) (from 0° to 50°).
- Padding with upholstery in many colours.
- Powder coated, very stable steel frame.
- Paper towel hanger (option).

## Technical data:

Number of sections:	2
Head rest angle adjustment [°]:	0/+50
Dimensions (L x W) [cm]:	193 x 61
Height [cm]	65
Max. paper width [cm]:	60
Weight [kg]:	16
Max. safe load [kg]:	150

## Accessories:



Paper roll hanger



Box dimen.: **60 x 59 x 24 cm** Weight: **+1,5 kg**

Stool for therapist **Flex (saddle stool)**: lightweight and comfortable stool for physicians and therapists with chrome-plated base and column.

## Features:

- Upholstery in many colours.
- Height adjustment.
- Tilting seat (front/back).
- Chrome-plated base.
- Rubberised castors.

## Technical data:

Height [cm]:	54-74
Tilting adjustment [°]:	15
Max. load [kg]:	150
Weight [kg]:	7,5
Base ø [cm]:	60
Seat dimension [cm]:	43x34

## Option:

- version with self breaking castors.



Box dimen.: **60 x 59 x 18 cm** Weight: **+1 kg**

Stool for therapist **Standard (round stool)**: lightweight and comfortable stool for physicians and therapists.

## Features:

- Upholstery in many colours.
- Height adjustment.
- Chrome-plated base.
- Rubberised castors.

## Technical data:

Height [cm]:	52-72
Max. load [kg]:	150
Weight [kg]:	6,5
Base ø [cm]:	60
Seat dimension ø [cm]:	36

## Option:

- version with self breaking castors.



# Qubiq, Medcart, Trolmed

Medical equipment trolleys

Excellent solution for electro- and physical therapy rooms. The construction of the trolley is designed specifically to accommodate the various types of medical equipments and allows easy access to them by appropriate arrangement of shelves in the space of the trolley.

The wheels ensure the mobility of the trolley and effortless maneuverability, and the brakes ensure stability and safety during treatments.

## Medical equipment trolley Qubiq

### Features:

- Compact shape.
- Stable.
- Made of aluminium, powder coated.
- 1 top shelf.
- 1 practical container for accessories.
- 2 drawers.
- Noiseless 4 transport wheels with brake.

### Technical data:

Dimensions (LxWxH) [cm]:	46 x 35 x 76
Max. load (max. safe trolley load with accessories) [kg]:	30
Drawer max. load [kg]:	5
Weight [kg]:	16,5



## Medical equipment trolley Medcart

### Features:

- Compact shape.
- Stable.
- Made of aluminium, powder coated.
- 2 shelves.
- Noiseless 4 transport wheels with brake.

### Technical data:

Dimensions (LxWxH) [cm]:	39 x 35 x 76
Max. load (max. safe trolley load with accessories) [kg]:	15
Shelf max. load [kg]:	6
Weight [kg]:	7



## Trolley for medical devices Trolmed type SPA-2

### Features:

- Compact shape.
- Powder coated steel.
- 2 shelves.
- 4 transport wheels; with brake.
- Removable 22l basket (optional accessory).
- Removable 5l basket (optional accessory).

optional accessory:



### Technical data:

Dimensions (WxDxH) [cm]:	54 x 43 x 89,5
Max. load (per shelf) [kg]:	25
Max. load (per basket / per small basket) [kg]:	5 / 2
Weight [kg]:	14





# Other equipment

As in this catalogue we presented only a part of our varied offer we encourage you to browse our full offer online:

[www.en.meden.com.pl](http://www.en.meden.com.pl)



You can find there:

- physiotherapy devices,
- devices for neurological rehabilitation
- many accessories as bolsters, wedges, stabilization belts, traction supports
- devices for supporting spine rehabilitation,
- hydrotherapy devices,
- patient lifts,
- neurological walkers,
- ... and many others.

# How to order easily

If our catalogue is of interest to you and you want to order any of presented products, please don't hesitate to contact us:



*European markets*

**Kamilla Modzelewska**  
Vice President

km@meden.com.pl  
+48 94 344 90 50  
+48 602 668 002



*CIS region, Ukraine, Georgia*

**Oksana Ardan**  
Key Account Specialist


oardan@meden.com.pl  
+48 94 347 10 67  
+48 607 566 302



*North/South America, Africa, Ireland, GB*

**Tomasz Górecki**  
Key Account Manager


tgorecki@meden.com.pl  
+48 601 646 866



*Australia, Asia and Middle East*

**Mariusz Rubik**  
Key Account Specialist  
& Training Consultant

mrubik@meden.com.pl  
+48 690 266 543



**Mateusz Baranowski**  
Training Consultant

mbaranowski@meden.com.pl  
+48 663 152 888



# Index

## A

Accessories 23

## C

ChiroSpace 51

## E

Elemental Ice 20

Elemental Fire 21

## F

Fiji P5/P7 53

## H

Hawaii S2/S4 53

## K

Kangoo 54

## L

LZM 57

## M

Medcart 58

## N

Nexus 48

## S

Safari Buffalo E2.C4-H 56

Safari Elephant B-S1 49

Safari Elephant B-S2 50

Safari Hydraulic Tables 46

Safari Jaguar P3 43

Safari Jaguar P5 44

Safari Leopard S3 41

Safari Leopard S5 42

Safari Puma E-S2 40

Safari Puma S2 38

Safari Puma S4 39

Safari Puma/Leopard (with CPR) 55

Stools 57

## T

Terapeuta Prestige E-S2 14

Terapeuta Prestige M-P3 8

Terapeuta Prestige M-P5 11

Terapeuta Prestige M-P7 13

Terapeuta Prestige M-S2 6

Terapeuta Prestige M-S3 7

Terapeuta Prestige M-S4 9

Terapeuta Prestige M-S5 10

Terapeuta Prestige M-S7 12

Terapeuta Prestige O-P5 16

Terapeuta Prestige O-S5 15

Theraspace E-S2 26

Theraspace M-S2 27

Theraspace M-S4 28

Theraspace Z3 29

Theraspace Z5 30

Theraspace Z7 31

Theraspace X4 32

Theraspace X6 33

Theraspace X8 34

Trolmed 58

## V

Vertimo Classic/Hi-Lo 52

## W

WaveMotion Ellipse 47

## Q

Qubiq 58



**Meden-Inmed Sp. z o.o.**

Wenedów 2, 75-847 Koszalin, Poland

Phone: +48 94 344 90 50, +48 94 347 10 40

Fax: +48 94 347 10 41

[international@meden.com.pl](mailto:international@meden.com.pl)

[www.en.meden.com.pl](http://www.en.meden.com.pl)

*We export to more  
than 100 countries:*



Follow us:

