

**Meden-Inmed Sp. z o. o.**

ul. Wenedów 2, 75-847 Koszalin, Poland

Phone.: +48 94 344 90 50, +48 94 347 10 40

Fax: +48 94 347 10 41

Mobile: +48 602 668 002

E-mail: [international@meden.com.pl](mailto:international@meden.com.pl)

[www.en.meden.com.pl](http://www.en.meden.com.pl)



**Theraspace**  
OMEGA UPHOLSTERY



Management  
System  
ISO 13485:2016

[www.tuv.com](http://www.tuv.com)  
ID 0000073249

For the production of our tables from the Theraspace line we use durable polyurethane upholstery.

Technical specification:

Surface: 100% Polyurethane  
Backing: 100% Polyester



Resistance to hydrolisis  
(8 years)



Does not contain chemical  
solvents



High scratch resistance



Respects the environment




Durability and performance  
>500 000 cycles (intensive usage)



Cleaning:

Any stain, spills or soiling should be cleaned up promptly with a neutral soap solution using a cloth in smooth circular movements.



OM006 WHITE



OM240 BEIGE



OM315 SKY



OM221 BLUE



OM314 CEMENT



OM076 OLIVE GREY



OM173 DARK BROWN



OM225 BLACK



OM469 LIMA



OM157 FOREST



OM053 RED



OM217 BORDEAUX